

VivoTrax™

Super-paramagnetic iron oxide tracer for pre-clinical Magnetic Resonance and Magnetic Particle Imaging

P/N MIVT01 5 vials

P/N MIVT01X 10 vials

Components: Each vial contains 1.1 mL of VivoTrax superparamagnetic iron oxide (SPIO) dextran magnetite nanoparticle solution




Format: Each vial contains ~5.5 mg/mL total iron in sterilized PBS
Endotoxins < 0.5 EU/mg Fe

Storage: Refrigerate at 2–8°C and away from sunlight

Properties: Hydrophilic colloidal solution of super-paramagnetic iron oxide coated with carboxydextran, with a mean hydrodynamic diameter of 62 nm

Appearance: Red-brown liquid

Stability: Guaranteed 6 months after shipment



A Molecular structure of VivoTrax tracer.

B Prussian blue staining of VivoTrax-labeled mesenchymal stem cells.³

C Tracking I.V. injected mesenchymal stem cells.³

For laboratory and animal research use only.

Warning: Not for human or animal therapeutic or diagnostic use. Make sure to comply with all laws and regulations governing research on animals.

Applications

Cell Tracking

Loading procedure, time of incubation and total uptake may differ based on cell line used. Please consider past publications for optimized protocols. The standard protocol below was developed based on bone marrow-derived human mesenchymal stem cell tracking³.

- Add 100 μL of VivoTrax tracer to 5 mL growth medium in T25 culture flasks (or at 100 $\mu\text{g}/\text{mL}$) containing adherent cells at 50% to 80% confluence to achieve desired Fe concentration.
- Incubate for 24 hours at 37°C with 5% CO_2 .
- After labeling, wash cultures twice with PBS to remove remaining free and surface-bound VivoTrax particles. Trypsinize cells, centrifuge and resuspend in PBS.
- Do not keep on ice for more than 2 hours before injection.

Blood Pool Imaging (4)*

Liver/Spleen	8 $\mu\text{mol Fe/kg}$	Bolus Injection
Perfusion	4–16 $\mu\text{mol Fe/kg}$	Bolus Injection
Angiography	10 $\mu\text{mol Fe/kg}$	Bolus Injection

MRI and Magnetic Particle Imaging

Have been shown to work in the majority of field strength MRI and MPI systems⁴, however please refer to manufacturer's imaging protocol.

VivoTrax is optimized for T2 and T2* weighted imaging MRI.

Relaxivities ($1 \text{ mM}^{-1} \text{ s}^{-1}$)¹

r_1 : 7.2 ± 0.1 (1.5 T and 37°C)

r_2 : 82.0 ± 6.2 (1.5 T and 37°C)

References

1. Capacity of human monocytes to phagocytose approved iron oxide MR contrast agents *in vitro*. Metz S, Bonaterra G, Rudelius M, Settles M, Rummeny E, Daldrup-Link H. *Eur Radiol*, 2004.
2. Impact of surface coating and particle size on the uptake of small and ultrasmall superparamagnetic iron oxide nanoparticles by macrophages. Saito S, Tsugeno M, Koto D, Mori Y, Yoshioka Y, Nohara S, Murase K. *Int J Nanomedicine*, 2012.
3. Quantitative Magnetic Particle Imaging Monitors the Transplantation, Biodistribution, and Clearance of Stem Cells *In Vivo*. Zheng B, von See MP, Yu E, Gunel B, Lu K, Vazin T, Schaffer DV, Goodwill PW, Conolly SM. *Theranostics* 2016; 6(3):291–301. doi:10.7150/thno.13728.
4. Superparamagnetic iron oxide contrast agents: physicochemical characteristics and applications in MR imaging. Wang YX, Hussain SM, Krestin GP. *Eur Radiol* 11:2319–2331.

For laboratory and animal research use only.

Warning: Not for human or animal therapeutic or diagnostic use. Make sure to comply with all laws and regulations governing research on animals.



980 Atlantic Avenue
Suite 102
Alameda, CA 94501
USA

Email: info@maneticinsight.com
magneticinsight.com

FOR RESEARCH USE ONLY
© Copyright 2016, Magnetic Insight. All rights reserved. MOMENTUM and VivoTrax are trademarks of Magnetic Insight.